



Polyclonal Anti-Integrin β2 (Sepharose Bead Conjugate)

Catalogue No. PA11224-S	Immunogen
Lot No. 08J01	A synthetic peptide mapping at the N-terminal of human Integrin β 2, different from the related mouse sequence by five amino acids.
Ig type: rabbit IgG	
Size: 100µg/vial	Purification Immunogen affinity purified.
Specificity Human, mouse, rat,. No cross reactivity with other proteins.	Formulation

50% slurry in

Recommended application

(Immunoprecipitation(IP)

PBS pH 7.2 with 0.01mg NaN₃a₃ preservative.

Storage Store at 4°C for frequent use.

Description:

This Antagene antibody is immobilized via covalent binding of primary amino groups to N-hydroxysuccinimide (NHS)-activated sepharose beads. It is useful for immunoprecipitation assays

BACKGROUND

The beta-2 integrin chain gene is designated ITGB2 and the leukocyte antigen has been designated CD18. The 3 alpha integrin chains associated individually with the beta-2 chain as a heterodimer have gene designations of ITGAL, ITGAM, and ITGAX, and leukocyte antigen designations of CD11A, CD11B, and CD11C, respectively. The expression of CD18 was increased in lymphoblastoid cells from persons with Down syndrome, consistent with the location of the gene on chromosome 211. The gene spans approximately 40 kb and contains 16 exons and all exon/intron boundaries conform to the GT/AG splicing consensus2. Furthermore,ITGB2 was constitutively clustered. Although it was expressed on the cell surface at normal levels and was capable of function following extracellular stimulation, it could not be activated via the "inside-out" signaling pathways3.

REFERENCE

1.Taylor, G. M.; Williams, A.; D'Souza, S. W.; Fergusson, W. D.; Donnai, D.; Fennell, J.; Harris, R. : The expression of CD18 is increased on trisomy 21 (Down syndrome) lymphoblastoid cells. Clin. Exp. Immun. 71: 324-328, 1988. 2.Weitzman, J. B.; Wells, C. E.; Wright, A. H.; Clark, P. A.; Law, S. K. A. : The gene organisation of the human beta-2 integrin subunit (CD18). FEBS Lett. 294: 97-103, 1991.

3.McDowall, A.; Inwald, D.; Leitinger, B.; Jones, A.; Liesner, R.; Klein, N.; Hogg, N. : A novel form of integrin dysfunction involving beta-1, beta-2, and beta-3 integrins. J. Clin. Invest. 111: 51-60, 2003.