



Product Information Sheet

Polyclonal Anti- LXR α (Magnetic Bead Conjugate)

Catalogue No. PA1246-M	Immunogen
Lot No. 09G01	A synthetic peptide corresponding to a sequence at the C-terminal of human LXR α , identical to the related rat and mouse sequence.
Ig type: rabbit IgG1	Purification
Size: 100 μ g/Vial	Immunogen affinity purified
Specificity	Contents
Rat.	Each vial contains 1mg/ml Magnetic Bead in PBS, pH 7.2, 0.05mg NaN ₃ .
No cross reactivity with other proteins.	Storage
	Store at 4°C for frequent use.
Recommended application	Description:
<i>Immunoprecipitation(IP)</i>	This Antagene antibody is immobilized by the covalent reaction of hydrazinonicotinamide-modified antibody with formylbenzamide-modified magnetic beads. It is useful for immunoprecipitation

BACKGROUND

LXRA is a tissue-specific cofactor that permits RXRA to function as a potent 9cRA receptor with a distinct target gene specificity.¹ It specifically interacts with RXRA in vivo to form a functional heterodimer in which RXRA is the ligand-binding subunit. Additionally, LXR activity is critical for physiologic lipid metabolism and transport.² LXRs are endogenous inhibitors of atherogenesis and are targets for therapeutic intervention in cardiovascular disease. Furthermore, LXRs and their ligands are negative regulators of macrophage inflammatory gene expression.³ LXR is also found that as a transcriptional switch that integrates hepatic glucose metabolism and fatty acid synthesis.⁴ The LXR-IDOL-LDLR axis defines a complementary pathway to sterol response element-binding proteins for sterol regulation of cholesterol uptake.⁵

REFERENCE

1. Willy, P. J.; Umesono, K.; Ong, E. S.; Evans, R. M.; Heyman, R. A.; Mangelsdorf, D. J. : LXR, a nuclear receptor that defines a distinct retinoid response pathway. *Genes Dev.* 9: 1033-1045, 1995.
2. Tangirala, R. K.; Bischoff, E. D.; Joseph, S. B.; Wagner, B. L.; Walczak, R.; Laffitte, B. A.; Daige, C. L.; Thomas, D.; Heyman, R. A.; Mangelsdorf, D. J.; Wang, X.; Lusis, A. J.; Tontonoz, P.; Schulman, I. G. : Identification of macrophage liver X receptors as inhibitors of atherosclerosis. *Proc. Nat. Acad. Sci.* 99: 11896-11901, 2002.
3. Joseph, S. B.; Castrillo, A.; Laffitte, B. A.; Mangelsdorf, D. J.; Tontonoz, P. : Reciprocal regulation of inflammation and lipid metabolism by liver X receptors. *Nature Med.* 9: 213-219, 2003.
4. Mitro, N.; Mak, P. A.; Vargas, L.; Godio, C.; Hampton, E.; Molteni, V.; Kreusch, A.; Saez, E. : The nuclear receptor LXR is a glucose sensor. *Nature* 445: 219-223, 2007.
5. Zelcer, N.; Hong, C.; Boyadjian, R.; Tontonoz, P. : LXR regulates cholesterol uptake through Idol-dependent ubiquitination of the LDL receptor. *Science* 325: 100-104, 2009.

For Research Use Only not for diagnostic and clinical use

Contact: Antagene, Inc. | Tel: 1 (866) 964-2589 | Fax: 1 (888) 225-1868 | Email: Info@antageneinc.com